

Otovent Adult



Category

Ear Care

Manufacturer

KESTREL MEDICAL LTD

Pack size

1 nose piece and 10 balloons

RRP

£14.99

Legal Status

MD



Uses

Drug-free treatment for Eustachian tube dysfunction.



Format and Ingredients

Nose piece



Dosage

Users inflate the balloon using their nostril while holding the other nostril closed. Use with each nostril up to 3 times a day. Blowing the balloon up creates enough pressure to open the Eustachian tube, making it easier for fluid to drain from the middle ear. Use up to 3 times daily until all fluid has drained away.