



Counter Intelligence PLUS | 2025

The training guide for pharmacy assistants



Canespro Fungal Nail Treatment Set

Uses

Removes the infected parts of the nail in 2-3 weeks.

Dosage

Step 1: Apply a thin layer of the ointment to the infected nail.

Step 2: Cover the treated nail with the plaster provided and keep covered for 24 hours.

Step 3: Remove the plaster and soak the nail in warm water for approximately 10 minutes. Gently remove the softened parts of the nail using the plastic scraper.

Repeat steps 1, 2 and 3 daily for 2-3 weeks.

Format and Ingredients

Treatment set containing:

10g tube Canespro Urea 40% ointment, 22 waterproof plasters, 1 plastic nail scraper,

Category:
Footcare

Manufacturer:
BAYER PLC

Pack size: 10g

RRP: £33.69

Legal Status: MD