

Counter Intelligence PLUS | 2025



The training guide for pharmacy assistants



Activa British Standard Class 1 Compression Hosiery

Uses

Provides 14-17mmHg of consistent graduated compression. Aids venous return and the prevention of venous disease progression. Effective for early varicose veins and in the immediate care of a wound on the lower limb in the absence of contraindications and appropriate onward referral and recommendations made.

Dosage

To be worn daily.

Format and Ingredients

Compression hosiery

Category:

Haemorrhoids and Circulatory

Manufacturer:

L&R MEDICAL UK LTD

Pack size: 1 pair

RRP: £15.90

Legal Status: MD





