



## Canespro Fungal Nail Treatment Set

### Uses

Removes the infected parts of the nail in 2-3 weeks. As well as Canespro Urea 40% ointment, the set includes 22 single use waterproof plasters and 1 plastic nail scraper.

### Dosage

Step 1: Apply a thin layer of the ointment to the infected nail.  
Step 2: Cover the treated nail with one of the plasters provided and keep covered for 24 hours.  
Step 3: Remove the plaster and soak the nail in warm water for approximately 10 minutes. Gently remove the softened parts of the nail using the plastic scraper.  
Repeat steps 1, 2 and 3 daily for 2-3 weeks.

### Format and Ingredients

Ointment containing:  
Urea 40%, Cera Alba, Lanolin, Petrolatum, Silica,

**Category:**  
Footcare, Skin Care and Antiseptics

**Manufacturer:**  
BAYER PLC

**Pack size:** Set

**RRP:** £38.29

**Legal Status:** MD