



DENTEK DENTAL FIRST AID KIT

Uses

Temporary dental repair kit to repair lost fillings or secure loose caps, crowns or inlays for up to 48 hours, until a dentist is seen. No mixing required.

Dosage

Follow step-by-step instructions on the pack. A dentist must be visited within 48 hours of use.

Format and Ingredients

Dental repair kit containing:

Anaesthetic / disinfectant Eugenol, Sesame Oil, Temporary filling material Zinc Oxide, Calcium Sulfate, Paraffinum Liquidum, Potassium Alum, Aluminium Sulfate, Aluminium Phosphate, Menthol, Eugenol, Iron Oxide Yellow,

Category:

Mouth Care

Manufacturer:

PRESTIGE CONSUMER HEALTH

Pack size:

RRP: £8.99

Legal Status: MD